Use of Xpert MTB/RIF and Xpert MTB/RIF Ultra on GeneXpert 10-colour instruments

WHO POLICY STATEMENT



2021



WHO POLICY STATEMENT

Use of Xpert MTB/RIF and Xpert MTB/RIF Ultra on GeneXpert 10-colour instruments

Use of Xpert MTB/RIF and Xpert MTB/RIF Ultra on GeneXpert 10-colour instruments: WHO policy statement

ISBN 978-92-4-004009-0 (electronic version) ISBN 978-92-4-004010-6 (print version)

© World Health Organization 2021

Some rights reserved. This work is available under the Creative Commons Attribution-NonCommercial-ShareAlike 3.0 IGO licence (CC BY-NC-SA 3.0 IGO; <u>https://creativecommons.org/licenses/by-nc-sa/3.0/igo</u>).

Under the terms of this licence, you may copy, redistribute and adapt the work for non-commercial purposes, provided the work is appropriately cited, as indicated below. In any use of this work, there should be no suggestion that WHO endorses any specific organization, products or services. The use of the WHO logo is not permitted. If you adapt the work, then you must license your work under the same or equivalent Creative Commons licence. If you create a translation of this work, you should add the following disclaimer along with the suggested citation: "This translation was not created by the World Health Organization (WHO). WHO is not responsible for the content or accuracy of this translation. The original English edition shall be the binding and authentic edition".

Any mediation relating to disputes arising under the licence shall be conducted in accordance with the mediation rules of the World Intellectual Property Organization (<u>http://www.wipo.int/amc/en/mediation/rules/</u>).

Suggested citation. Use of Xpert MTB/RIF and Xpert MTB/RIF Ultra on GeneXpert 10-colour instruments: WHO policy statement. Geneva: World Health Organization; 2021. Licence: <u>CC BY-NC-SA 3.0 IGO</u>.

Cataloguing-in-Publication (CIP) data. CIP data are available at http://apps.who.int/iris.

Sales, rights and licensing. To purchase WHO publications, see <u>http://apps.who.int/bookorders</u>. To submit requests for commercial use and queries on rights and licensing, see <u>https://www.who.int/copyright</u>.

Third-party materials. If you wish to reuse material from this work that is attributed to a third party, such as tables, figures or images, it is your responsibility to determine whether permission is needed for that reuse and to obtain permission from the copyright holder. The risk of claims resulting from infringement of any third-party-owned component in the work rests solely with the user.

General disclaimers. The designations employed and the presentation of the material in this publication do not imply the expression of any opinion whatsoever on the part of WHO concerning the legal status of any country, territory, city or area or of its authorities, or concerning the delimitation of its frontiers or boundaries. Dotted and dashed lines on maps represent approximate border lines for which there may not yet be full agreement.

The mention of specific companies or of certain manufacturers' products does not imply that they are endorsed or recommended by WHO in preference to others of a similar nature that are not mentioned. Errors and omissions excepted, the names of proprietary products are distinguished by initial capital letters.

All reasonable precautions have been taken by WHO to verify the information contained in this publication. However, the published material is being distributed without warranty of any kind, either expressed or implied. The responsibility for the interpretation and use of the material lies with the reader. In no event shall WHO be liable for damages arising from its use.

Contents

Acknowledgementsiv
Abbreviations and acronymsv
WHO policy statementvi
Background1
GeneXpert 10-colour versus 6-colour instruments1
Summary of methods2
Summary of results
TAG meeting outcome4
Remarks4
Implementation considerations5
Further research5
References
Annex 1: List of participants7
Technical Advisory Group (TAG) members (in alphabetical order)7
Annex 2: Declaration of interests10
Annex 3: Meeting agenda15
Web Annex: Study report

Acknowledgements

The development of the policy statements was led by Nazir Ahmed Ismail, Alexei Korobitsyn and Carl-Michael Nathanson, with support from Matteo Zignol, and under the overall direction of Tereza Kasaeva, Director of the World Health Organization (WHO) Global Tuberculosis (TB) Programme (WHO/GTB). The WHO/GTB gratefully acknowledges the support and contributions of the following individuals:

Technical Advisory Group

Heidi Albert, Foundation for Innovative New Diagnostics (FIND), the global diagnostics alliance, South Africa; Khalide Azam, Southern Africa TB Health System Support Project, East, Central and Southern Africa Health Community, United Republic of Tanzania; Daniela Cirillo, San Raffaele TB Supranational Reference Laboratory (SRL), Milan, Italy; Christopher Coulter, Queensland Mycobacterium Reference Laboratory (SRL) and Communicable Diseases Branch, Queensland Health, Brisbane, Australia; Claudia Maria Denkinger, Division of Tropical Medicine, Heidelberg University Hospital, Germany; Patricia Hall, Centers for Disease Control and Prevention (CDC), Atlanta, United States of America; Nguyen Van Hung, Department of Microbiology, National Tuberculosis Reference Laboratory, Viet Nam; Farzana Ismail, Centre for Tuberculosis, National Institute for Communicable Diseases (NICD)/ National Health Laboratory Service (NHLS); SRL, Johannesburg, South Africa; Irina Lyadova, Laboratory of Cellular and Molecular Basis of Histogenesis, Koltzov Institute of Developmental Biology of Russian Academy of Sciences, Russian Federation; Florian Maurer, National Reference Centre for Mycobacteria, Research Center Borstel, Germany; Sandeep Meharwal, FHI360, Thailand; Vithal Prasad Myneedu, South Asian Association for Regional Cooperation (SAARC) TB and HIV/AIDS Centre, Nepal; Mark Nicol, University of Western Australia, Perth, Australia; Alaine Umubyeyi Nyaruhirira, Management Sciences for Health, South Africa; Madhukar Pai, McGill International TB Centre, McGill University, Canada; Paulo Redner, National Reference Laboratory for Tuberculosis, Oswaldo Cruz Foundation, Brazil; Sadia Shakoor, Departments of Pathology and Pediatrics, Aga Khan University Hospital Karachi, Pakistan; Siva Kumar Shanmugam, Department of Bacteriology, National Institute for Research in Tuberculosis, Indian Council of Medical Research, India; Xin Shen, Division of Tuberculosis and HIV/AIDS Prevention, Shanghai Municipal Center for Disease Control and Prevention, China; Thomas Shinnick, Independent Consultant, Atlanta, GA, United States of America; Sabira Tahseen, National Tuberculosis Control Programme, Ministry of National Health Services, Regulations and Coordination, Government of Pakistan, Islamabad, Pakistan; Yanlin Zhao, National Tuberculosis Control and Prevention Center, Chinese Centers for Disease Control and Prevention, China.

Technical review team

Elisa Tagliani, San Raffaele SRL, Milan, Italy; **Shaheed Vally Omar**, Centre for Tuberculosis, National Institute for Communicable Diseases (NICD)/ National Health Laboratory Service (NHLS); SRL, Johannesburg, South Africa;

Observers

Draurio Cravo Neto Barreira, Unitaid, Geneva, Switzerland; **Fatim Cham-Jallow**, The Global Fund to Fight AIDS, Tuberculosis and Malaria, Geneva, Switzerland; **Brian Kaiser**, Global Drug Facility, Stop TB Partnership, Geneva, Switzerland; **Kaiser Shen**, United States Agency for International Development (USAID), United States of America; **Andrei Mosneaga**, Stop TB Partnership, Geneva, Switzerland.

Funding

This product was developed with support from USAID and the Russian Federation.

Abbreviations and acronyms

AIDS	acquired immunodeficiency syndrome
Ct	cycle threshold
DR-TB	drug-resistant tuberculosis
FIND	Foundation for Innovative New Diagnostics
GTB	Global TB Programme
GXP	GeneXpert instrument
GXP6	GeneXpert 6-colour instrument
GXP10	GeneXpert 10-colour instrument
HIV	human immunodeficiency virus
МТВ	Mycobacterium tuberculosis
МТВС	Mycobacterium tuberculosis complex
PQ	prequalification
RIF	rifampicin
SRL	WHO TB Supranational Reference Laboratory
TAG	Technical Advisory Group
ТВ	tuberculosis
Tm	melting temperature
Ultra	Xpert MTB/RIF Ultra
USA	United States of America
WHO	World Health Organization
XDR	extensively drug-resistant

WHO policy statement

Introducing improved, rapid and more accurate diagnostic tools¹ is critical to achieving the global targets towards ending the tuberculosis (TB) epidemic. The World Health Organization (WHO) recommended the use of Xpert[®] MTB/RIF and Xpert MTB/RIF Ultra (Ultra) (Cepheid, Sunnyvale, United States of America [USA]) for the detection of TB and rifampicin-resistant TB in 2010 and 2017, respectively.² Both tests are widely implemented as initial tests for patients with presumptive TB and are performed on GeneXpert instruments with 6-colour optics. In 2021, WHO recommended the class of low complexity automated nucleic acid amplification tests (NAATs) to detect resistance to amikacin, ethionamide, fluoroquinolones and isoniazid.³ The first-in-class test is the Xpert MTB/XDR (Cepheid, Sunnyvale, USA). In contrast to the Xpert MTB/RIF and Ultra, this test requires an instrument with 10-colour optics and cannot be performed on the existing 6-colour instrument systems. The performance of Xpert MTB/RIF and Ultra on the new GeneXpert 10-colour instruments has not been independently assessed. Having a single instrument that could be used to detect TB and resistance to first-line and second-line drugs would simplify workflow and facilitate implementation.

To evaluate the evidence on the performance of Xpert MTB/RIF and Xpert Ultra on GeneXpert 10colour instruments, WHO convened a meeting of the Technical Advisory Group (TAG) on Tuberculosis Diagnostics and Laboratory Strengthening on 5–6 October 2021.⁴ This document provides background information and describes the available evidence and discussions by the TAG.

Following review of the evidence and advice from the TAG, WHO makes the following policy statements:

- 1. Xpert MTB/RIF and Ultra cartridge performance on the GeneXpert 10-colour instrument is comparable to that of the GeneXpert 6-colour instrument for detection of TB and rifampicin resistance.
- 2. Current WHO recommendations for Xpert MTB/RIF and Ultra cartridge use on GeneXpert 6colour instruments are also valid for their use on GeneXpert 10-colour instruments.

The guidance provided should facilitate procurement and uptake of these technologies and improve patient care. The statements above should be read in the context of the remarks and implementation considerations detailed in this report. In addition, further research questions are proposed that seek to address data gaps and inform models to improve effective implementation of the tests. The current WHO recommendations on the use of the Xpert MTB/RIF and Ultra on the GeneXpert 6-colour instrument and the use of the Xpert MTB/XDR on the GeneXpert 10-colour instrument are unchanged and remain valid. All products recommended by WHO are automatically eligible to be included in the

预览已结束, 完整报告链接和二维码如下:



https://www.yunbaogao.cn/report/index/report?reportId=5 23467